

BIPOLAR COAGULATION FORCEPS

angled on shanks; blunt; Ø 1.0mm

17.0cm (6.8"); flat plug

SKU: 022-107-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [curved](#)