

BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

16.0cm (6.4"); flat plug



SKU: 022-110-16

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.67.05 |
| SPECIFICATIONS | sharp , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.67.05 |
| SPECIFICATIONS | sharp , straight |