

## BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

18.0cm (7"); flat plug



**SKU:** 022-110-18

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [sharp](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [sharp](#), [straight](#)