

## BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

22.0cm (8.8"); flat plug



**SKU:** 022-110-22

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.23.03
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.23.03
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>