

BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

24.0cm (9.6"); flat plug



SKU: 022-110-24

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.24.03
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.24.03
SPECIFICATIONS	sharp , straight