

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 1.0mm

16.0cm (6.4"); flat plug



SKU: 022-112-16

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.67.06
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.67.06
SPECIFICATIONS	blunt , straight