

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 1.0mm

20.0cm (8"); flat plug



SKU: 022-112-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.70.06 |
| SPECIFICATIONS | blunt , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.70.06 |
| SPECIFICATIONS | blunt , straight |