

## BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 1.0mm

22.0cm (8.8"); flat plug



**SKU:** 022-112-22

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Aesculap</b>	GK970
<b>Medicon</b>	88.23.12
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Aesculap</b>	GK970
<b>Medicon</b>	88.23.12
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">straight</a>