

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 1.0mm

22.0cm (8.8"); flat plug



SKU: 022-112-22

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Aesculap	GK970
Medicon	88.23.12
SPECIFICATIONS	blunt , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Aesculap	GK970
Medicon	88.23.12
SPECIFICATIONS	blunt , straight