

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug



SKU: 022-112-24

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.24.12
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.24.12
SPECIFICATIONS	blunt , straight