

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

20.0cm (8"); flat plug



SKU: 022-114-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.16
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.16
SPECIFICATIONS	blunt , straight