

## BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

22.0cm (8.8"); flat plug



**SKU:** 022-114-22

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**SPECIFICATIONS** [blunt](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**SPECIFICATIONS** [blunt](#), [straight](#)