

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

22.0cm (8.8"); flat plug



SKU: 022-114-22

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)