

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

24.0cm (9.6"); flat plug



SKU: 022-114-24

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [blunt](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [blunt](#), [straight](#)