

BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.3mm

24.0cm (9.6"); flat plug



SKU: 022-100-24

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [sharp](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [sharp](#), [straight](#)