

BIPOLAR COAGULATION FORCEPS

bajonett; angled up.; pointed; Ø 0.3mm

22.0cm (8.8"); flat plug



SKU: 022-121-22

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [sharp](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [sharp](#)