

## BIPOLAR COAGULATION FORCEPS

bajonett; angled up.; pointed; Ø 0.3mm

22.0cm (8.8"); flat plug



**SKU:** 022-121-22

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [sharp](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [sharp](#)