

BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 1.0mm

18.0cm (7"); flat plug



SKU: 022-123-18

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)