

BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 1.0mm

20.0cm (8"); flat plug



SKU: 022-123-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.38
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.38
SPECIFICATIONS	angled , blunt