

BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug



SKU: 022-123-24

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)