

## BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 2.0mm

18.0cm (7"); flat plug



**SKU:** 022-125-18

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)