

BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.3mm

15.0cm (6"); flat plug

SKU: 022-101-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.07
SPECIFICATIONS	angled , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.07
SPECIFICATIONS	angled , sharp