

## BIPOLAR COAGULATION FORCEPS

bajonett; angled down; pointed; Ø 0.3mm

24.0cm (9.6"); flat plug

**SKU:** 022-131-24

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [sharp](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [sharp](#)