

## BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 1.0mm

20.0cm (8"); flat plug



**SKU:** 022-133-20

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                           |
| <b>TYPE</b>           | <a href="#">Bipolar Coagulation Forceps</a>    |
| <b>Medicon</b>        | 88.70.39                                       |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">blunt</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                           |
| <b>TYPE</b>           | <a href="#">Bipolar Coagulation Forceps</a>    |
| <b>Medicon</b>        | 88.70.39                                       |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">blunt</a> |