

BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 1.0mm

22.0cm (8.8"); flat plug



SKU: 022-133-22

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)