

## BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug



**SKU:** 022-133-24

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)