

BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.3mm

20.0cm (8"); flat plug



SKU: 022-101-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
Medicon 88.20.07
SPECIFICATIONS [angled](#), [sharp](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
Medicon 88.20.07
SPECIFICATIONS [angled](#), [sharp](#)