

BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.3mm

15.0cm (6"); flat plug



SKU: 022-100-15

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.05
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.05
SPECIFICATIONS	sharp , straight