

## BLOOD PRESSURE MANOMETERS

wall mounted; calibrated

with velcro cuff



**SKU:** 100-580-00

**Categories:** [Diagnostics](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Blood Pressure Manometers</a>
<b>KLS Martin</b>	17-994-00

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Blood Pressure Manometers</a>
<b>KLS Martin</b>	17-994-00