

BLOUNT KNEE RETRACTOR

45mm; 2 teeth

26.0cm (10.4")

SKU: 190-576-26

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
KLS Martin	23-565-26
AUTHORS	BLOUNT
Medicon	77.12.52



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
KLS Martin	23-565-26
AUTHORS	BLOUNT
Medicon	77.12.52