

BLOUNT KNEE RETRACTOR

38mm; blunt

26.0cm (10.4")

SKU: 190-575-26

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
KLS Martin	23-563-26
AUTHORS	BLOUNT
Medicon	77.12.51



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
KLS Martin	23-563-26
AUTHORS	BLOUNT
Medicon	77.12.51