

BÖHLER BONE RONGEUR

straight; double convertible

15.0cm (6")


SKU: 196-307-01

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Rongeurs
KLS Martin	23-876-15
AUTHORS	BOEHLER
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Rongeurs
KLS Martin	23-876-15
AUTHORS	BOEHLER
SPECIFICATIONS	straight