

## BONE WIRE

Ø 0.4mm; coil of 10m



**SKU:** 190-795-04

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Twisting Instruments</a>
<b>Aesculap</b>	KY004
<b>KLS Martin</b>	22-527-04
<b>Medicon</b>	76.01.04

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Twisting Instruments</a>
<b>Aesculap</b>	KY004
<b>KLS Martin</b>	22-527-04
<b>Medicon</b>	76.01.04