

BONNEY DRESSING FORCEPS

straight; serrated

18.0cm (7.2"); width of jaw 3.2 +/-0.2

SKU: 020-120-18

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD078 |
| AUTHORS | WAUGH |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD078 |
| AUTHORS | WAUGH |
| SPECIFICATIONS | serrated , straight |