

## BORCHARDT FLAT TREPHINE

with eye; shaft 5mm

Ø 3.5mm



**SKU:** 190-284-35

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-194-35
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-194-35
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>