

## BORCHARDT FLAT TREPHINE

with eye; shaft 5mm

Ø 5.0mm



**SKU:** 190-284-50

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-194-50
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-194-50
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>