

BORCHARDT TWIST DRILL

short working end; shaft 5mm

Ø 2.5mm



SKU: 190-283-25

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-192-25 |
| AUTHORS | BORCHARDT |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-192-25 |
| AUTHORS | BORCHARDT |