

BORCHARDT TWIST DRILL

short working end; shaft 5mm

Ø 2.5mm



SKU: 190-283-25

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-25
AUTHORS	BORCHARDT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-25
AUTHORS	BORCHARDT