

## BORCHARDT TWIST DRILL

short working end; shaft 5mm

Ø 2.8mm



**SKU:** 190-283-28

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-192-28
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-192-28
<b>AUTHORS</b>	<a href="#">BORCHARDT</a>