

BORCHARDT TWIST DRILL

short working end; shaft 5mm

Ø 3.5mm



SKU: 190-283-35

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-35
AUTHORS	BORCHARDT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-35
AUTHORS	BORCHARDT