

BORCHARDT TWIST DRILL

long working end; shaft 5mm

Ø 3.2mm



SKU: 190-282-32

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-190-32
AUTHORS	BORCHARDT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-190-32
AUTHORS	BORCHARDT