

BORCHARDT TWIST DRILL

long working end; shaft 5mm

Ø 6.0mm



SKU: 190-282-60

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-190-60 |
| AUTHORS | BORCHARDT |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-190-60 |
| AUTHORS | BORCHARDT |