

## BORCHARDT TWIST DRILL

long working end; shaft 5mm

Ø 2.0mm



**SKU:** 190-282-20

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |                             |
|-------------------|-----------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>        |
| <b>TYPE</b>       | <a href="#">Hand Drills</a> |
| <b>KLS Martin</b> | 22-190-20                   |
| <b>AUTHORS</b>    | <a href="#">BORCHARDT</a>   |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                             |
|-------------------|-----------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>        |
| <b>TYPE</b>       | <a href="#">Hand Drills</a> |
| <b>KLS Martin</b> | 22-190-20                   |
| <b>AUTHORS</b>    | <a href="#">BORCHARDT</a>   |