

## BOZEMANN-DOUGLAS DRESSING FORCEPS

s-shaped; serrated; big Ring

26.0cm (10.4")

**SKU:** 030-327-26

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	EO085
<b>KLS Martin</b>	14-315-26
<b>AUTHORS</b>	<a href="#">BOZEMANN</a> , <a href="#">DOUGLAS</a>
<b>Medicon</b>	19.82.26
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	EO085
<b>KLS Martin</b>	14-315-26

<b>AUTHORS</b>	<a href="#">BOZEMANN</a> , <a href="#">DOUGLAS</a>
<b>Medicon</b>	19.82.26
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>