

BOZEMANN DRESSING FORCEPS

straight; serrated

26.0cm (10.4")




SKU: 030-320-26

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	EO080
AUTHORS	BOZEMANN
Medicon	19.80.26
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	EO080
AUTHORS	BOZEMANN
Medicon	19.80.26

SPECIFICATIONS

[serrated](#), [straight](#)