

## BOZEMANN DRESSING FORCEPS

straight; serrated

26.0cm (10.4")

**SKU:** 030-320-26

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	EO080
<b>AUTHORS</b>	<a href="#">BOZEMANN</a>
<b>Medicon</b>	19.80.26
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	EO080
<b>AUTHORS</b>	<a href="#">BOZEMANN</a>
<b>Medicon</b>	19.80.26

## SPECIFICATIONS

[serrated](#), [straight](#)