

## BRAND TENDON HOLDING FORCEPS

straight; 2x3 teeth

19.0cm (7.6")



**SKU:** 195-241-19

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF326
<b>KLS Martin</b>	23-676-19
<b>AUTHORS</b>	<a href="#">BRAND</a>
<b>Medicon</b>	78.10.19
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF326
<b>KLS Martin</b>	23-676-19

## AUTHORS

Medicon

## SPECIFICATIONS

[BRAND](#)

78.10.19

[straight](#), [toothed](#)