

BRAND TENDON HOLDING FORCEPS

straight; 2x3 teeth

24.0cm (9.6")



SKU: 195-241-24

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF328
KLS Martin	23-676-24
AUTHORS	BRAND
Medicon	78.10.25
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF328
KLS Martin	23-676-24

AUTHORS

Medicon

SPECIFICATIONS

[BRAND](#)

78.10.25

[straight](#), [toothed](#)