

BRAND TENDON HOLDING FORCEPS

straight; 2x3 teeth

15.0cm (6")



SKU: 195-241-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF324
KLS Martin	23-676-15
AUTHORS	BRAND
Medicon	78.10.15
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF324
KLS Martin	23-676-15

AUTHORS

Medicon

SPECIFICATIONS

[BRAND](#)

78.10.15

[straight](#), [toothed](#)