

BROWN NEEDLE HOLDER

straight; diagonally serrated

18.0cm (7.2")



SKU: 050-096-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
KLS Martin	20-344-18
AUTHORS	BROWN
Medicon	10.15.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
KLS Martin	20-344-18
AUTHORS	BROWN
Medicon	10.15.18

SPECIFICATIONS

[serrated](#), [straight](#)